

Description

OT Bioelettronica offers a double-layer adhesive medical sheet, for EMG acquisitions, to be interposed between the HD matrices and intact skin.

The sheet is present in 3 different configurations with different IED interelectrode distance. The reduced thickness and adhesiveness guarantee optimal positioning and contact.

CND Z12100485

Materials

Health

Technical details

Medical layer

Tape	Double coated PE 0,16mm
Carrier	Trasparent PE 0,08mm
Adhesive	Acrylate designed for medical use

Industrial layer

Tape	PE
Adhesive	Acrylate

Packaging

Envelope	Paper ATICELCA 501
----------	--------------------

Biocompatibility (UNI EN ISO 10993)	Si
--	----

Latex	No
-------	----

Phthalates	No
------------	----

Generic code	KITXXMMXXYY
--------------	-------------

Sheet size	340x250mm
------------	-----------

Hole size	3mm
-----------	-----


IED distance	5mm, 8mm, 10mm
--------------	----------------

Rows x Cols number	67x49
--------------------	-------

	42x31
--	-------

	34x25
--	-------

How to use the double adhesive foams



Visit www.otbioelettronica.it for video tutorials